

VISUPAC Star System requirements

VISUPAC Server	
Operating System	Windows 2000, Windows 2003 Server, Windows XP
CPU	Pentium III or higher
RAM	1 GB
Network	100 Mbit; TCP/IP
Storage	Hard drive 300 GB + NAS 2 TB for archiving (Structure and size of the storage system may be adapted to the respective usage conditions.)
VISUPAC Viewer	
Operating System	Windows 2000, Windows XP Pro, Windows Vista Business & Ultimate
CPU	Pentium IV 2 GHz or higher
RAM	2 GB (768 MB minimal)
Graphic Card	32 Bit color, minimal resolution 1280 x 1024
Network	100 Mbit, TCP/IP
Storage	Hard drive >20 GB
VISUPAC Star connection interfaces	
Transfer images/reports to VISUPAC	<ul style="list-style-type: none"> • DICOM file import • Generic file import (Jpg, Tiff, BMP, GIF, Doc, Xml, PDF) • Print2Visupac (printer driver for Win 2000 & Win XP based devices)
Transfer patient master data from VISUPAC to connected instruments	<ul style="list-style-type: none"> • DICOM modality worklist scheduler
Transfer patient master data from *PMS and *EMR systems to VISUPAC	<ul style="list-style-type: none"> • SMP Basic (proprietary text file based patient data import) • DICOM modality worklist interface

*PMS: patient management system *EMR: electronic medical record (system)

Publication No: 000000-1696-586
The contents of the brochure may differ from the current status of approval of the product in your country. Please contact our regional representative for more information.
Subject to change in design and scope of delivery and as a result of ongoing technical development. Printed on elemental chlorine-free bleached paper. PUBLICIS I/2009.
© 2008 by Carl Zeiss Meditec AG. All copyrights reserved.

VISUPAC Star Storage and Review System

Know it all.



Carl Zeiss Meditec AG
Goeschwitzer Str. 51-52
07745 Jena
GERMANY

Phone: +49 36 41 22 03 33
Fax: +49 36 41 22 01 12
info@meditec.zeiss.com
www.meditec.zeiss.com

Carl Zeiss Meditec Inc.
5160 Hacienda Drive
Dublin, CA 94568
USA

Phone: +1 925 557 41 00
Fax: +1 925 557 41 01
info@meditec.zeiss.com
www.meditec.zeiss.com



VISUPAC Star™ Storage and Review System

Your connection to enhanced workflow, better patient care.

Knowledge is power: the power to provide your patients with the best possible care in the most efficient way.

With the VISUPAC Star system, you can know virtually everything there is to know about a given patient's diagnostic test results with your Carl Zeiss Meditec instruments, all with a level of immediacy and convenience that saves time and lets you keep your priority on providing the highest level of patient care.

Connectivity is the Key.

The VISUPAC Star system offers centralized storage for test reports and clinical images generated by your diagnostic instruments. It also provides easy access to the images and reports from secure and easy-to-use workstations in your exam lanes.

The VISUPAC Star system enables you to arrive at even more confident conclusions with greater convenience and more effective use of existing resources through:

Central Viewing

The VISUPAC Star system will now provide you with the ability to view all diagnostic images and reports without having to search for outputs on the individual modalities.

Central Management

The VISUPAC Star system will store all patient diagnostic images and reports on a central server for ease of data management and backup.

Minimizing Charting Errors

The new patient scheduler eliminates the need for typing patient information into each diagnostic instrument, thereby saving time, ensuring patient data is consistent, and minimizing charting errors.

With the VISUPAC Star system, you can count on:

- The ability to capture and store diagnostic images and/or reports from most Carl Zeiss Meditec instruments and from many 3rd party instruments.
- Immediate access to diagnostic images and reports in all exam lanes.
- The capability to set up different user accounts and to track image and report modifications by specific users
- Interfaces with DICOM enabled patient management and electronic medical record systems.

For more information on how to connect with greater confidence and convenience, contact your Carl Zeiss Meditec representative today and we will be happy to consult with you and your IT staff to determine how best to integrate VISUPAC Star into your everyday clinical workflow.

