

IOL		Diopter Range & Injector / Cartridge
4-Haptic MICS	CT SPHERIS® 209M	0.0 to +32.0 D
	CT ASPHINA® 409M	AT.Shooter A1-2000 injector / VISCOJECT™-BIO 1.8 Cartridge Set
	CT ASPHINA 509M	or VISCOJECT-BIO 1.8 Injector Set
	AT LISA® 809M	or A6 injector & *AT.Smart cartridge
		or AT Shooter A2-2000 injector & *AT.Smart cartridge
	AT TORBI® 709M	-10.0 to +28.0 D sphere; +1.0 to +4.0 D cylinder
		AT.Shooter A1-2000 injector / VISCOJECT™-BIO 1.8 Cartridge Set
		or VISCOJECT-BIO 1.8 Injector Set
		or A6 injector & *AT.Smart cartridge
		or AT Shooter A2-2000 injector & *AT.Smart cartridge
		-10.0 to +28.0 D sphere; +4.5 to +12.0 D cylinder
		VISCOJECT-BIO 2.2 Injector Set
		or AT.Shooter A1-2000 injector / VISCOJECT-BIO 2.2 cartridge Set
		+28.5 to +32.0 D sphere; +1.0 to +12.0 D cylinder
		VISCOJECT-BIO 2.2 Injector Set
	or AT.Shooter A1-2000 injector / VISCOJECT-BIO 2.2 cartridge Set	
AT LISA toric 909M	-10.0 to +24.0 D sphere; +1.0 to +12.0 D cylinder	
	AT.Shooter A1-2000 injector / VISCOJECT™-BIO 1.8 Cartridge Set	
	or VISCOJECT-BIO 1.8 Injector Set	
	or A6 injector & *AT.Smart cartridge	
	or AT Shooter A2-2000 injector & *AT.Smart cartridge	
	+24.5 to +32.0 D sphere; +1.0 to +12.0 D cylinder	
	VISCOJECT-BIO 2.2 Injector Set	
	or AT.Shooter A1-2000 injector / VISCOJECT-BIO 2.2 cartridge Set	
AT LISA tri toric 939M	-10.0 to +28.0 D sphere; +4.5 to +12.0 D cylinder and	
	+28.5 to +32.0 D sphere; +1.0 to +12.0 D cylinder	
	VISCOJECT-BIO 2.2 Injector Set	
	or AT.Shooter A1-2000 injector / VISCOJECT-BIO 2.2 cartridge Set	
AT LARA toric 939M	-4.0 to +32.0 D SE; +4.5 to +12.0 D cylinder	
	VISCOJECT-BIO 2.2 Injector Set	
	or AT.Shooter A1-2000 injector / VISCOJECT-BIO 2.2 cartridge Set	
<b>BLUEMIXS® Preloaded IOLs</b>		
CT ASPHINA 409MP	0.0 to +32.0 D	BLUEMIXS 180 injector
CT ASPHINA 509MP		
AT LISA 809MP		
AT LISA tri 839MP		
AT LARA 829MP		
AT TORBI 709MP	-10.0 to +28.0 D sphere; +1.0 to +4.0 D cylinder	BLUEMIXS 180 injector
AT LISA toric 909MP	-10.0 to +24.0 D sphere; +1.0 to +4.0 D cylinder	BLUEMIXS 180 injector
AT LISA tri toric 939MP	-10.0 to +28.0 D sphere; +1.0 to +4.0 D cylinder	BLUEMIXS 180 injector
AT LARA toric 929MP	-8.0D to +32.0D spherical equivalent; +1.0 to +4.0 D cylinder	

		BLUEMIXS 180 injector
4-Haptic	CT SPHERIS® 204	<b>-10.0 to +30.0 D</b> AT.Shooter A1-2000 / VISCOJECT-BIO 2.2 Cartridge Set <i>or</i> VISCOJECT-BIO 2.2 Injector Set
	CT ASPHINA® 404	<b>+31.0 to +45.0 D (CT SPHERIS 204)</b> <b>+31.0 to +42.0 D (CT ASPHINA 404)</b> ACCUJECT™ 3.0 injector set
	CT SPHERIS 203	<b>+8.0 to +32.0 D</b> VISCOJECT-BIO 2.2 injector set
	<b>SkyJet® Preloaded IOLs</b>	
	CT SPHERIS 203P	<b>+8.0 to +26.0 D</b> SkyJet 2.8 injector <b>+27.0 to +30.0 D</b> SkyJet 3.2 injector
	CT ASPHINA 603P	<b>+10.0 to +30.0 D</b> SkyJet 2.8 injector
1-Pc C-Loop	CT 47S CT 37A AT LISA 801	<b>0.0 to +30.0 D</b> AT.Shooter A1-2000 / VISCOJECT-BIO 2.2 Cartridge Set <i>or</i> VISCOJECT-BIO 2.2 Injector Set <b>+31.0 to +40.0 D</b> ACCUJECT™ 3.0 injector set
	CT 47LC	<b>-10.0 to +30.0 D</b> AT.Shooter A1-2000 / VISCOJECT-BIO 2.2 Cartridge Set <i>or</i> VISCOJECT-BIO 2.2 Injector Set <b>+31.0 to +40.0 D</b> ACCUJECT™ 3.0 injector set
	<b>Preloaded hydrophobic IOLs</b>	
	CT LUCIA 601P/PY	<b>+0.0 to +24.0 D</b> ACCUJECT 2.0 injector set
	CT LUCIA 201P	<b>+24.5 to +30.0 D</b> ACCUJECT 2.2 injector set
	CT LUCIA 611P/PY	<b>+0.0 to +24.0 D</b> ACCUJECT 2.0 injector set <b>+24.5 to +30.0 D</b> ACCUJECT 2.2 injector set <b>+30.5 to +34.0 D</b> ACCUJECT 2.4 injector set
	CT LUCIA 211P	
	3 Piece	CT LUCIA 202
CT LUCIA 602 <sup>4</sup>		<b>+4.0 to +34.0 D</b> R28 injector / Z28 cartridge
	CT 13A	<b>+10.0 to +30.0 D</b> Forceps
	CT 27SF	<b>+5.0 to +30.0 D</b> Forceps

## Injector and Cartridge Type and Incision Size

Injector	Cartridge	Injector Type	Incision Size
AT.Shooter A2-2000	AT Smart Cartridge <sup>1</sup>	Reusable	1.8mm
AT.Shooter A1-2000	VISCOJECT-BIO 1.8 Cartridge set <sup>1,2</sup>	Reusable	
VISCOJECT-BIO 1.8 Injector Set <sup>2</sup>		Disposable	
Single Use Injector A6	*AT.Smart Cartridge Set <sup>1</sup>	Disposable <sup>2</sup>	
BLUEMIXS 180		Disposable, Pre-Loaded	
AT.Shooter A1-2000	VISCOJECT-BIO 2.2 Cartridge set <sup>1,2</sup>	Reusable	≥2.2mm
VISCOJECT-BIO 2.2 Injector Set <sup>2</sup>		Disposable	
ACCUJECT 2.0		Disposable	2.2mm
ACCUJECT 2.2		Disposable	2.4mm
SkyJet 2.8		Disposable, Pre-Loaded	≥2.8 mm ≤ 3.0mm
OD 655 DualTec kit 3-P Acrylic 2.8		Disposable	
R28 injector <sup>4</sup>	Z 28 cartridge <sup>4</sup>	Reusable	
ACCUJECT 3.0		Disposable	

<sup>1</sup> Soft silicone cushion is included in the cartridge set

<sup>2</sup> VISCOJECT - BIO or "classic" versions

<sup>3</sup> VISCOJECT 1.8 ≥ 1.8mm with docking technique

<sup>4</sup> Only available in the USA